

```
In [1]: %matplotlib inline
```

```
In [2]: # uNMR COMS is the first code operating the HMC832 PLL on the uNMR Board
```

```
# Set-Up
# 1. Connect the FTDI USB-RS485-WE to the uNMR board
# 2. Connect 5.2 - 5.5Vdc to the uNMR board
# 3. Programme the TMS320F28335PTPQ with the required code
# 4. Click "Cell" above and then click "Run All"
# 5. Measure PLL output initial value is
```

```
In [3]: #%matplotlib inline                                     # Polts plots in scrip

import time                                                    # Timer units seconds
import math                                                    # math functions
import numpy as np                                             # multi-dimensional arrays/
    matrices
import pandas as pd                                           # building block for data a
    nalysis
import minimalmodbus                                           # modbus RTU serial comms m
    odule
import matplotlib.pyplot as plt                               # Python plotting library

instrument = minimalmodbus.Instrument('com3', 1)              # port name, slave address
    (in decimal)
instrument.serial.baudrate = 57600                             # Baud
instrument.serial.bytesize = 8
instrument.serial.parity   = 'E'                               # Parity

instrument.serial.stopbits = 1
instrument.serial.timeout  = 0.05                              # seconds

instrument.mode = minimalmodbus.MODE_RTU                      # rtu or ascii mode
```

```
In [4]: vtune = instrument.read_registers(6,1)                 # read case 6
        print (vtune)
```

```
[3074]
```

```

In [7]: # Temperature sensor, board SPI, routine
# This routine requires "import numpy as np" for np.linspace
# This routine requires "import time"

temp_read = [] # array temp_read

count_max = 10 # sample count max
count_min = 0 # sample count min
yaxis_min = 25 # yaxis min
yaxis_max = 32 # yaxis max
print 'Count Temp (degC)' # Column header

for i in range(count_max):
    print '{:02f}'.format(i) # Prints every result
    time.sleep(.2) # Loop delay 0.2s
    val = instrument.read_register (7, 1) # read modbus reg 7 decima
    l 1
    val2 = val * 0.3125 # raw data to deg C conver
    sion
    print ' ', i, ' ', val2 # prints column header
    temp_read += [val2]


plt.plot( temp_read ) # Plot results
plt.title('Temperature vs Time')
plt.xlabel('Sample Count')
plt.ylabel('TEMP (deg C)')
plt.grid( True )

plt.axis([count_min, count_max, yaxis_min, yaxis_max])
plt.show()

#print temp_read

```

Count	Temp (degC)
0.000000	
0	29.5
1.000000	
1	29.5
2.000000	
2	29.5
3.000000	
3	29.5
4.000000	
4	29.5
5.000000	
5	29.46875
6.000000	
6	29.5
7.000000	
7	29.5
8.000000	
8	29.46875
9.000000	
9	29.46875



In [ ]: